# SHIMADZU MANUAL BY JJ

JAN.JAROS2@VUT.CZ

# OBSAH

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TRAPEZIUMX-V - NiTi\_cylinder\_d8mm.vmax

File Edit View Test Tools Hardware Window Help

System Sensor **Testing** Specimen Data Processing Chart Report

Copy Insert Delete Clear ← →

	Area1	Area2	Area3	Area4	Area5	Area6	Area7	Area8	Area9	Area10
Act.	Down	Up	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF
	Force	Force								
	1000.000	1000.000								
Change point	Channel	Channel								
	Force	Force								
	40000	100								
	N	N								
GetData	Stress Disp. Strain	Stress Disp. Strain								
Samplings	200msec	Same as prev. area								
Loop	→ 50Cycle	←								

Disp. Origin: Start  
 Auto control starting point: 0.1 %/FS  
 Initial: 1 mm/min

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TRAPEZIUMX-V - NiTi\_cylinder\_d8mm.vmax

File Edit View Test Tools Hardware Window Help

System Sensor **Testing** Specimen Data Processing Chart Report

Copy Insert Delete Clear ← →

	Area1	Area2	Area3	Area4	Area5	Area6	Area7	Area8	Area9	Area10
Act.	Down	Up	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF
	Force	Force								
	1000.000	1000.000								
Change point	Channel	Channel								
	Force	Force								
	40000	100								
GetData	Stress Disp. Strain	Stress Disp. Strain								
	200msec	Same as prev. area								
Samplings										
Loop	→ 50Cycle	←								

Disp. Origin: Start  
 Auto control starting point: 0.1 %/FS  
 Initial: 1 mm/min

End Setti...  
 Break Detection:  
 Sensitivity: 10.0 %  
 Level/%FS: 0.02 %  
 Level/%MAX: 50.0 %  
 Break and Limit Action:  Stop  
 Break Detection start point: 0.035 %

Pre-Test

< Back Next >

Save a Method file  
 Test with this method  
 Finish  
 Cancel

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TRAPEZIUMX-V - NiTi\_cylinder\_d8mm.vmax

File Edit View Test Tools Hardware Window Help

System Sensor **Testing** Specimen Data Processing Chart Report

Copy Insert Delete Clear ← →

	Area1	Area2	Area3	Area4	Area5	Area6	Area7	Area8	Area9	Area10
Act.	Down	Up	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF
	Force	Force								
	1000.000	1000.000								
Change point	Channel	Channel								
	Force	Force								
	40000	100								
GetData	Stress Disp. Strain	Stress Disp. Strain								
	200msec	Same as prev. area								
Loop	→ 50Cycle	←								

Disp. Origin: Start  
 Auto control starting point: 0.1 %/FS  
 Initial: 1 mm/min

End Setti...  
 Pre-Test

Break Detection:  
 Sensitivity: 10.0 %  
 Level/%FS: 0.02 %  
 Level/%MAX: 50.0 %

Break and Limit Action:  
 Stop

Break Detection start point: 0.035 %

< Back Next >


Save a Method file  
 Test with this method  
 Finish  
 Cancel

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TRAPEZIUMX-V - NiTi\_cylinder\_d8mm.vmx

File Edit View Test Tools Hardware Window Help

System Sensor Testing Specimen **Data Processing** Chart Report

Normal Mode

Setup D.P. Items.  
Settings of the data processing options.

- 1 Click on the formula icons
- 2 Set parameters
- 3 Set pass/fail options
- 4 Select the Statistics
- 5 Define any custom formulas

Save a Method file

Test with this method

Finish

Cancel

< Back Next >

**Control - Compression**

**Data Processing**

Data Processing	Print
GetData1_Force_Lp1	<input checked="" type="checkbox"/>
GetData1_Stress_Lp1	<input checked="" type="checkbox"/>
GetData1_Displ_Lp1	<input checked="" type="checkbox"/>
GetData1_Strain_Lp1	<input checked="" type="checkbox"/>
GetData1_Time_Lp1	<input checked="" type="checkbox"/>
GetData2_Force_Lp1	<input checked="" type="checkbox"/>
GetData2_Stress_Lp1	<input checked="" type="checkbox"/>
GetData2_Displ_Lp1	<input checked="" type="checkbox"/>
GetData2_Strain_Lp1	<input checked="" type="checkbox"/>
GetData2_Time_Lp1	<input checked="" type="checkbox"/>
GetData1_Force_Lp2	<input checked="" type="checkbox"/>
GetData1_Stress_Lp2	<input checked="" type="checkbox"/>
GetData1_Displ_Lp2	<input checked="" type="checkbox"/>
GetData1_Strain_Lp2	<input checked="" type="checkbox"/>
GetData1_Time_Lp2	<input checked="" type="checkbox"/>

Edit Delete

**Statistics**

- Average
- Standard Deviation
- Maximum
- Minimum
- Range
- Median
- Variation
- 3Sigma
- Average+ 6 Sigma
- Average- 6 Sigma

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
TRAPEZIUMX-V - AZ61\_WAAM\_plocha\_3.vmav

File Edit View Test Tools Hardware Window Help

System Sensor Testing Specimen Data Processing Chart Report

Preview Open Template Save Template Page Setup

ReportDesigner

Tensile test\_AI3D

Key Word	Product Name		
Test File Name	Method File Name		
Report Date	15.02.2024	Test Date	AZ61_WAAM_plocha_3.vmav
Test Mode	Single	Test Type	Tensile
Speed	0.00021/sec	Shape	Plate
No of Batches:	2	Qty/Batch:	5

Name	Elastic	YS1_Force	YS1Rp	Max_Rm	Break_Strength	Max_Fm	Break_Force	Break_Ext1(Strain)
Parameters	Ext. 1(Strain) 0.005 - 0.025 %	0.2 %	0.2 %	Calc. at Entire Areas	Sensitivity : 10	Calc. at Entire Areas	Sensitivity : 10	Sensitivity : 10
Unit	GPa	N	N/mm2	N/mm2	N/mm2	N	N	%

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Test with this method

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# TESTOVÁNÍ

The screenshot shows a software interface with a 'Test Wizard' dialog box. The 'Specimen' tab is selected, displaying the following information:

- Material: Metal
- Shape: Rod
- No of Batches: 3
- Qty/Batch: 2
- Size Unit: mm

Below the specimen information, there are two tables for 'Cfm Of Caliper Input'.

Name	Diameter [D]	Height [H]	Data/Constant		
			Type	Diameter	
1- 1	38	8.0000	10.0000	tube	2
1- 2	38	8.0000	10.0000	tube	2
2- 1	56	8.0000	10.0000	tube	2
2- 2	56	8.0000	10.0000	tube	2
3- 1	76	8.0000	10.0000	tube	2
3- 2	76	8.0000	10.0000	tube	2

The background shows a graph of Stress (MPa) vs Ext. 1 (Strain(%)) and a 'Chart1' window.

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**Test Wizard - Specimen**

Material: Metal | No of Batches: 3 | Size Unit: mm  
 Shape: Rod | Qty/Batch: 2

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